**Mountain View
Medical Supply**
5076 W. 58th Ave.
Arvada, CO 80002
303-455-1300
800-873-7121
MVMSinc.com

VH SonicMist Ultrasonic Nebulizer System

Model 11-520

**READ THIS INSTRUCTION MANUAL
CAREFULLY BEFORE USE**

INSTRUCTION MANUAL • ENGLISH & ESPAÑOL

STOP!

**PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING
COMPONENTS BEFORE USING YOUR NEBULIZER**

UltraSonic Unit

Mouthpiece #1

Nasal Piece

Child and Adult
Masks

Mouthpiece #2

Extension Tube

Air Filter

5 Medication Cups

Unit Base

Battery Pack


Auto Adapter

AC Adapter

Storage Bag

Instruction Manual
FAQ Information Sheet
Replacement Parts
Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR
INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF
PURCHASE. CONTACT CUSTOMER CARE AT 866y326y1313.**

TABLE OF CONTENTS

Introduction	4
Care and Safety Information	5-6
Product Features	7
Nebulizer Mask/Mouthpiece Selection	8
Power Option Use and Selection	9-10
Using Your Nebulizer	11-14
Cleaning and Disinfection	15-16
Changing the Filter	17
Troubleshooting	18
Warnings and Cautions	19
Product Specifications	20
Warranty	21
Instrucciones en Español.....	23-43

**Toll-Free Customer Care Help Line: 1-866-326-1313
Monday - Friday 8:30 a.m. - 4:30 p.m. CST**

Distributed for:
Veridian Healthcare
1175 Lakeside Drive Gurnee, IL
60031
www.veridianhealthcare.com
Made in China
#93-1069 07/12
©2012 Veridian Healthcare, LLC

INTRODUCTION

Thank you for purchasing this Veridian Healthcare VH SonicMist Ultrasonic Nebulizer. Please read this instruction manual completely before using your unit to ensure safe and effective use, as well as long life of the product.

Ultrasonic nebulizers are used to efficiently deliver medication in the form of a mist to the lungs. Prescribed medicine is placed in the nebulizer and the ultrasonic transponder converts the liquid medication into an aerosol mist which is then inhaled. This method of delivery is beneficial in the treatment of asthma and other respiratory ailments, as the medicine is administered directly to the bronchial passages of the lungs to provide quick and effective relief. With proper use and care, your new nebulizer will provide you with many years of safe, effective breathing therapy treatments.

Federal law allows nebulizers to be sold by a physician's prescription only. These devices should be used only at the direction of a doctor or healthcare professional and should be used only as intended and described in this manual.

CARE AND SAFETY INFORMATION

Improper care or use of your nebulizer may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the nebulizer's efficacy and long life.

General Cautions and Warnings

- Adult supervision is recommended for use with children
- To avoid electrical shock keep unit away from water, never immerse the adapter or unit in liquid; do not reach for a unit that has fallen into water—immediately unplug the unit
- Do not handle the unit or power cord(s) with wet hands
- Do not use the unit if it has any damaged parts (including power cord or plug), if it has been submersed in water or dropped
- The unit should not be used where flammable gas, oxygen or aerosol spray products are being used
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Do not leave the unit unattended when plugged in or operating
- Do not shake the unit when in operation
- Opening the cover during operation may cause injury
- Disconnect the unit from the electrical outlet before cleaning, filling and after each use
- Only use included attachments and accessories; do not use attachments from other nebulizer models
- Masks and mouthpieces should be replaced at least once a year
- Do not disassemble the unit or attempt to repair it

Operating Cautions and Warnings

- Do not allow nebulized medication to come in contact with eyes; nebulizer treatments are intended to be inhaled orally or nasally only
- Change the water after each use
- The maximum capacity of the medicine cup is 7.5 ml; do not overfill

CARE AND SAFETY INFORMATION

- Always follow physician instruction for volume of medication used per treatment and necessary treatment time
- Never use this unit while operating a vehicle
- If you experience any discomfort or an abnormality occurs, stop using the unit immediately
- The unit should not be used for more than 30 minutes at a time, after use wait 30 minutes before beginning a new treatment

Storage Cautions and Warnings

- Never store the unit in direct sunlight, high temperature or humidity
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use; it is also suggested that the battery pack be disconnected from the unit when not in use

Cleaning Cautions and Warnings

- Always check the selected mouthpiece, medication cup and air filter before each use; dirty or worn parts should be cleaned or replaced as necessary
- Never immerse the unit in water to clean as it may damage the unit
- Disconnect the unit from the electrical outlet or remove battery pack before cleaning
- Clean all necessary parts after each treatment as instructed in this guidebook; if the unit has not been used for an extended period, all parts should be cleaned prior to use

Battery Cautions and Warnings

- Included battery pack is suitable for use with VH SonicMist only, do not use with other units
- Do not use any other power source except those provided; unsupported power sources may damage unit
- Do not store the battery in direct sunlight, high temperature or high humidity
- Avoid placing the battery on metal or wet surfaces or inside a metal box; this may short-circuit the battery's positive and negative contact causing the battery to heat or explode

PRODUCT FEATURES

Illustrations of the nebulizer mechanics and detailed use instructions are included throughout this manual.

Air Filter Cover

Cover

POWER Button

Airflow Control/
Fan Speed Button

Battery Pack

NEBULIZER MASK/MOUTHPIECE SELECTION

This unit includes five styles of masks and mouthpieces, allowing for user preference in their desired medication delivery method. Replacement parts are available for order on the enclosed order form, by calling Customer Care at 866-326-1313 or at your local medical supply dealer.


Face mask (adult and child size included): face masks require use of the included extension tube. Select the appropriate size for the patient - child or adult - and place over the patients nose and mouth. Secure in place with the elastic strap around the patient's head and hold the nebulizer in place with a free hand.

Mouthpiece style #1: put the mouthpiece into the mouth, holding gently between the teeth and close the lips around the mouthpiece; hold the nebulizer in place with a free hand.


Mouthpiece style #2: this style requires use of the extension tube. Put the mouthpiece into the mouth, holding gently between the teeth and close the lips around the mouthpiece; hold the nebulizer in place with a free hand.

Nasal piece: gently insert the nasal pieces into the nose; hold the nebulizer in place with a free hand.

Proper use and assembly of each mouthpiece is further explained later in the 'USING YOUR NEBULIZER' section of this manual.

Adult Face Mask

Mouthpiece Style #1

Mouthpiece Style #2

Nasal Piece

POWER OPTIONS USE AND SELECTION


This unit includes three convenient power options:

Rechargeable
Battery Pack

AC Adapter

DC Auto
Adapter

Replacement power option parts are available for order on the enclosed order form, by calling Customer Care at 866-326-1313 or at your local medical supply dealer.


Note: To remove the nebulizer body from the battery pack, press the ribbed release button as shown and pull up on the nebulizer body.

Rechargeable Battery Pack:

Note: Charge for 24 hours before first use; charge for 4-6 hours after initial charge. Battery can be charged with AC adapter plug or DC auto adapter.

1. Remove the battery pack from the housing.
2. Insert the adapter into the battery pack power port and plug the plug end into a wall socket.

A green light will illuminate on the battery and remain on while the unit is charging; this light does not necessarily indicate that the battery is fully charged, always follow the recommended charging times as noted above.

It is recommended that the battery be removed after each use and when the nebulizer will not be used for an extended period of time.

POWER OPTIONS USE AND SELECTION

AC Adapter:

1. Be certain the battery pack is removed from the housing.
2. Insert the AC adapter directly into the underside of the housing, Fig 1. For hands-free operation, use with the included base, Fig 2.
3. Plug the AC adapter plug end into a wall socket.
4. Always unplug the adapter before cleaning and storage.


Fig. 1

DC Auto Adapter:

1. Be certain the battery pack is removed from the housing.
2. Insert the DC adapter directly into the underside of the housing, Fig 3.
3. Plug the adapter directly into the power source in the vehicle.
4. Always unplug the adapter before cleaning and storage.


Fig. 2

Fig. 3

Never leave the unit plugged in to a power source when not in use.

USING YOUR NEBULIZER

Before the first use of the nebulizer or if the unit has not been used for an extended period of time, it should be cleaned thoroughly as described in the **Cleaning and Disinfection** section of this manual.



The patient receiving treatment should be in a comfortable seated or slightly reclining position as treatments take several minutes.

1. Determine your power option and mouthpiece/mask style selection and gather your parts prior to starting the treatment.
2. Attach your selected power source.
3. Remove the clear nebulizer cover, Fig 1. Remove the medication cup if one is in place.
4. Fill the water reservoir with either tap or distilled water to the indicated level between the two lines within the reservoir, Fig 2. Use fresh water for every treatment.

Note: The recommended water temperature range is 68°F-104°F/20°C-40°C

The use of cool water (below 68°F/20°C) will result in a lower rate of nebulization.

The use of warm water (above 104°F/40°C) will result in a higher rate of nebulization.

USING YOUR NEBULIZER

5. Insert the medication cup as shown, Fig 1.
6. Pour the medication into the nebulizer cup following the physician prescribed dosage, Fig 2. The medication cup capacity is 7.5mL; overfilling the medication cup will reduce the nebulization rate.
Note: If your medication dosage volume is 2ml or less, or if the medication is very dense, the medication should be diluted with distilled water so that the volume is a minimum of 4 ml.
7. Replace the clear nebulizer cover, Fig 3.
8. Plug the power cord into the electrical outlet if necessary.


Fig. 1


Fig. 2






Fig. 3

Proper use of the mask or mouthpiece is important in helping to ensure effective and efficient medication delivery. Some pieces require an extension tube (included) for proper use. An extension tube is flexible and can be adjusted for a more comfortable position during treatment.

USING YOUR NEBULIZER


9. Attach the selected mouthpiece to the nebulizer and position the mouthpiece at the patient's mouth.

Remember that masks should cover the patient's mouth and nose; mouthpieces should be in the mouth with the lips closed around the end of the mouthpiece, Figs 4-7.

USING YOUR NEBULIZER

10. Press the POWER button  on the nebulizer. Nebulization begins at a slower rate and will become stronger.
11. To increase or decrease the nebulization rate, press the fan button  until the desired level is reached. The lights on the housing will illuminate in correlation to the nebulization rate: 3 lights = highest rate, 1 light = lowest rate.
12. The nebulizer should remain vertical during treatment for proper nebulization of the medication and to prevent spilling of the liquid.
13. Ultrasonic treatment times average six minutes. The viscosity of the medication and airflow setting may result in longer treatment times.
14. When the treatment is complete, turn the unit off and unplug it from the outlet.
15. Always dispose of any remaining medication.
16. Allow the unit to cool. Disassemble, clean and store your nebulizer. Be certain that the nebulizer pieces are completely dry before storage.

It is important to complete the entire prescribed dose of medication for maximum effectiveness. Consult with your prescribing physician or healthcare professional for information on your specific medication therapy treatment plan.

CLEANING AND DISINFECTING PROCEDURES

The nebulizer and mouthpiece should be thoroughly cleaned after each use unless a different cleaning regimen is specified by your healthcare professional; additional hygienic precautions may be required for use with high-risk patients. Never attempt a subsequent medication treatment without first cleaning your unit and mouthpiece.

Cleaning Procedures:

1. Ensure that the power cord is disconnected from the power outlet and that the battery pack is disconnected from the nebulizer housing body.
2. Remove the clear lid from the nebulizer body, Fig 1.
3. Wash the mouthpiece, fan, medicine cup and air filter cover in warm water with mild detergent. Rinse thoroughly, Fig 2.
4. Allow all pieces to air dry before use or storage. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.
5. Rinse the water reservoir with warm water, do not use detergent. As needed, wipe out the reservoir with a damp cloth. Do not replace the lid or medication cup until the unit is completely dry.
6. The body of the unit and the outside of the battery pack can be wiped down with a soft cloth; if necessary, use a damp soft cloth with a non-abrasive cleaner. Never immerse the nebulizer body or battery in liquid.


Fig. 1


Fig. 2

CLEANING AND DISINFECTING PROCEDURES

Disinfecting Procedures:

It is recommended that the parts be disinfected monthly if the unit is used on a regular medication therapy schedule.

1. Remove the clear lid from the nebulizer body.
2. Rinse the mouthpiece, fan, medicine cup and air filter cover in hot water. Mix one part white vinegar with three parts distilled water, in a volume sufficient to fully immerse the pieces.
3. Submerge pieces in the vinegar mixture for thirty minutes.
4. Rinse everything thoroughly with water.
5. Allow all pieces to air dry completely, or dry them thoroughly with a clean, soft towel. Do not store any parts wet as remaining moisture creates a risk of bacterial growth.

Never use soap pads or powder cleanser to clean any part of the unit. Never use any cleaning products that may be harmful when swallowed or inhaled.


It is recommended that mouthpieces and masks be replaced periodically depending on the frequency of therapy treatments or when damaged or cracked. To order replacement parts, use the included order form, call Customer Care at 866-326-1313 or contact your local dealer.

CHANGING THE FILTER

- Never operate the unit without a filter.
- Replace the filter after every ten uses, or sooner if it turns gray, appears contaminated or becomes wet.
- Use only VH SonicMist filters; additional filters are available for order using the included order form, calling Customer Care at 866-326-1313 or contacting your local dealer.
- Do not substitute cotton or any other materials for the filters.
- Do not attempt to clean or reuse a filter.

To replace the filter, follow these steps:

1. Change the filter only when the unit is off, unplugged and completely cool.
2. Gently pull the filter cap out of the cover; if necessary, use a non-metal tool to help lift up cap.
3. Remove the filter from the filter cap.
4. Place a new filter in the filter cap.
5. Return the filter cap to its place, ensuring that it is firmly in place.

TROUBLESHOOTING

The nebulizer produces no aerosol, or a small amount of aerosol.

- There may be too much or too little medicine in the nebulizer. The amount of liquid should be between 0.5 ml and 7.5 ml.
- The medication cup may be damaged; check and replace.
- The water reservoir may be dirty; wipe clean and attempt again.
- The nebulizer should be held vertically during use.
- The medication may not be suitable for nebulization. Check with your physician.
- The filter may be in need of replacement

Some medicine remains in the nebulizer at the end of my treatment.

This is normal. Never keep leftover medication.

Treatments with the mask takes longer than treatments with the mouthpiece.

This is normal. Less medication is delivered per breath with the mask because the aerosol is combined with ambient air through the holes in the mask.

Unit does not operate when the POWER button is pressed.

If using the AC adapter or auto adapter, check connections between the adapter and the unit and the adapter and the outlet. If using the battery, check that the battery is properly connected. Charge the battery.

Fan is not working properly.

Check that the fan is firmly pushed onto the motor shaft.

Fan speed light changes from red to green during use.

The battery needs to be charged. Nebulization will slow and then cease.

WARNING AND CAUTIONS

Nebulizer therapy is intended for use with physician prescribed medication only. Over the counter medication, essential/aromatherapy oils or alcohol are not suitable for inhaling using a nebulizer. Highly viscous substances may damage the unit and make it unsuitable for use. Never use nebulizer therapy for anything other than physician prescribed respiratory care therapy.

Use common sense electrical caution when operating the device. This nebulizer runs on DC power. Always turn off the power and unplug the device before cleaning or if the unit is not working properly.

PRODUCT SPECIFICATIONS

Name	VH SonicMist Ultrasonic Nebulizer
Model Number	11-520
Medication Capacity	Maximum 7.5 ml
Particle Size	0.5 to 10 μm
Ultrasonic Frequency	2.5MHz+200kHz (-100 kHz)
Average Nebulization Rate	$\geq 0.2 \text{ ml / minute}$
Operating Conditions	Temperature: 41°F - 104°F (5°C - 40°C); Humidity: 10% - 90%
Unit Dimensions (W x H x D)	3.35" x 2.35" x 5" (without battery) 3.35" x 2.35" x 6.25" (with battery)
Unit Weight	8 oz. (without battery) 13.8 oz. (with battery)
Power Connection	12V DC
Included Accessories	Rechargeable battery pack, 12-volt AC adapter, DC auto adapter, child mask, adult mask, extension tube, nasal piece, 2 mouthpiece styles, 5 medication cups, 1 air filter, base

Specifications are subject to change without notice.

3yYEAR LIMITED WARRANTY

Veridian Healthcare guarantees that its VH SonicMist Ultrasonic Nebulizer will be free from manufacturing defects under normal use for a period of three-years. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the nebulizer unit. This warranty is voided if the nebulizer unit is misused or abused in any manner.

If the nebulizer fails to operate within the warranty period while under the original purchaser owns the product, return it postage prepaid along with \$5.00 for return shipping and insurance to: Veridian Healthcare, Attn: Repair Department, 1175 Lakeside, Gurnee, IL, 60031; include your name, address, phone number and brief description of the problem with your nebulizer. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. Veridian Healthcare will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you.


This warranty constitutes the warrantor's only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. Veridian will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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
**Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

Distributed for:
Veridian Healthcare
1175 Lakeside Drive Gurnee, IL
60031
www.veridianhealthcare.com

Made in China
#93-1069 07/12
©2012 Veridian Healthcare, LLC

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