



USAGE GUIDELINES FOR SWIFT SLIDER



Mountain View Medical Supply
5076 W. 58th Ave.
Arvada, CO 80002
303-455-1300
800-873-7121
MVMSinc.com






BREATHABLE PATIENT REPOSITIONER



STEP 1


Place Swift Slider on bottom sheet with blue, low-friction taffeta panel down.

Note: Ensure brakes of bed are on.



STEP 2a

Tuck in flaps, once Swift Slider has been placed.



STEP 2b

Place incontinence pad on top of Swift Slider (if needed).



STEP 3

Position patient centrally on white Swift Slider surface. The central panel should support patient from the neck to just below the buttocks.

Note: Place pillow at head of bed to avoid accidental head injury.



STEP 4

Ensure there is a caregiver on each side of bed.

Note: When possible, bed height should be at caregivers' upper thigh level.




STEP 5

Un-tuck flaps.



STEP 6


Bend patient's knees to protect heels from being dragged across linens.



STEP 7

Hold Swift Slider with two hands - One at shoulder height and one at hip height.

To move patient, gather Swift Slider as close as possible to patient's body. Reposition by laterally shifting body weight.



STEP 7

Alternate repositioning method: Rabbit Ears technique



STEP 8

Once patient has been repositioned, smooth out Swift Slider and tuck flaps securely under mattress.



USAGE GUIDELINES FOR SWIFT SLIDER

Note: If the head of the bed needs to be elevated, raise the knee gatch to prevent the patient from sliding down the mattress.

Position the patient's hips in the bend of the bed. If hips go above that zone, patient will slide down by force of gravity.

Please remember:

We strongly advise against using Swift Slider for lifting.



Swift