

4-Layer Compression Bandage

Indications: Manage venous leg ulcers and related conditions. This system consists of four bandages that, when applied correctly, provide compression for seven days. Multi-layered system simultaneously absorbs, protects, and provides graduated sustained compression. Not made with natural rubber latex.

Contraindications: Determine the Ankle Brachial Pressure Index (ABI) using a Doppler ultrasound.

- An ABI less than 0.8 may indicate the presence of significant arterial disease.
- Do not use this Bandaging System if the patient's ABI is less than 0.8 or if arterial disease is suspected.
- While using this product, if the patient develops pain, numbness, color changes, cool toes or foot, the patient should seek immediate medical attention.
- Do not use this Bandaging System on patients with an ankle circumference less than 18 cm.
- Apply extra padding as needed over bony prominences of the ankle and/or tibia to avoid pressure necrosis.

Change Frequency: Change dressing as directed by a healthcare professional.

Recommended Secondary Dressings: n/a



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY	EXUDATE LEVEL*
Compression Bandage System, 4-Layer, Non-Sterile	1008	1/pk, 8pk/cs	N/A

EXUDATE LEVEL KEY: ◐◐◐ = Light ◐◐◐ = Moderate ◐◐◐ = Heavy

*Exudate level indicates highest level of exudate that the dressing can accommodate



Mountain View Medical Supply
 5076 W. 58th Ave.
 Arvada, CO 80002
 303-455-1300 or 800-873-7121
 MVMSinc.com

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.