

## REBOUND® AIR WALKER

The Rebound Air Walker features a host of details designed to provide added protection, superior compression, simplified application and usage, improved comfort and, ultimately, successful patient outcomes.

**Mountain View Medical Supply**  
**5076 W. 58th Ave.**  
**Arvada, CO 80002**  
**303-455-1300 or 800-873-7121**  
**MVMSinc.com**



### INDICATIONS FOR USE

#### REBOUND AIR WALKER

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase. If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.

#### REBOUND AIR WALKER LOW TOP

Foot and ankle conditions that may benefit from motion restriction, such as:

- Stable foot fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase.

### FEATURES AND BENEFITS

- Full circumferential shell with additional semi-soft toe bumper for frontal toe protection
- Soft Flex Edge™ reduces pressure points
- Full rocker bottom is especially engineered to be low and wide to help promote a natural gait
- Stretch-to-fit shell for optimal fit
- Shell can be trimmed for additional stretch in the calf
- Cushioning and breathable softgood provides comfort
- Linked medial/lateral bladders for easy inflation and equal compression
- Integrated pump and air release valve for easy access

### REQUIRED MEASUREMENTS

- Shoe size

