

STAND ASSIST PATIENT LIFT

Through assisted standing, the lift is an aid to rehabilitation and encourages strength and stability from the patient, keeping the patient active and engaged. The patient lift raises patient from sitting position to standing position, useful for the frequent daily transfer between bed, chair and bathroom. This lift is engineered to provide safety and comfort for both the resident and the caregiver.





Rechargeable Battery

Direct or standalone charging options are available.

Height Adjustable Leg Cushion

Patient Handles

Patient handles on both sides of the

boom to enhance stability of the patient.

The patient leg support can be easily adjusted to suit different patient body types. The soft ergonomic cushion provides additional comfort for the patient.



The double roller front castors improves the lift stability.



SWL: 180kg / 400lb

Adjustable Base

Extended shifter handle enables easy and quick opening and closing of the lift base.

Lockable Castors

The rear castors can be locked to keep the lift stationary when unused.

Anti-slip Foot Stand

Dedicated patient foot stand with anti-slip finishing.



Stand-up Sling



Transfer Sling (LPS23)



Charging Station (LPR10)

Height at Sling	Height at Sling	Base Height	Weight of Lift	Weight Capacity
Hook Up (Max)	Hook Up (Min)	(Furniture Clearance)		(Safe Working Load)
1778 mm / 70"	853 mm / 33.6"	114 mm / 4.5"	55 kg / 122lb	180 kg / 400lb



PATIENT SLINGS


Costcare provides an extensive range of slings to be paired with the Costcare patient lifts. The variety of slings offer solutions to all lifting circumstances.



Divided-leg Sling (LPS21)



Transfer Sling (LPS23)



Stand-up Sling (LPS24)



Full Body Sling (LPS26)



Full Body Mesh With Commode Sling (LPS28)

Sling	Sizes				
Model	X-Large	Large	Medium	Small	
LPS21	0	0	0	0	
LPS23		0	0		
LPS24	0	0			
LPS26	0	0	0	0	
LPS28		0	0		